



- Compact housing, color white RAL 9010
- Ceramic heating element with 600 W
- 4-step electronic temperature adjustment
- 15-minute alarm (audible and visible) with automatic power reduction to safety level of 20%.
- Built-in two-stage LED lighting
- Power failure alarm
- Applied standards: EN 60601-1, EN 60601-2-21 and EN 60601-1-2
- **With double-joint arm, ceiling anchorage and height adjustment. Dimensions as per drawing**
- Weight: 30.0 kg

More details:

Datasheet Ceratherm 600-3 Nr. 15.249.X_

